

SIMULATOR MODELS

Aorta | Carotis | Leg

Abdominal model for "open aortic surgery"

The VASCULAR INTERNATIONAL (VI) abdominal model works independently and can be used anywhere. The waterfilled aortic inlay, under pulsatile pressure, makes a lifelike training of the AAA repair possible. All arteries from the aorta are open and can be anastomosed.




Abdominal model prepared for open AAA repair



View of the abdominal aorta after tube graft interposition



Preparation for a bifurcated graft



Fixation of the aortic inlay below the bifurcation

VASCULAR INTERNATIONAL offers lifelike aortic inlays which fit perfectly into the VI abdominal model. All arteries from the aorta are open and can be anastomosed. The sewability of the inlays is very lifelike – for interrupted as well as continuous sutures. The density of the anastomosis and with it the operation result can easily be checked.



All in one equipment

Each simulator comes with a full pump set including pump control and cable-tie kit. It can easily be transported in a box.

SALES PRICE incl. pump set, 1 aorta, cable tie kit CHF **6.850,-** (ex works Switzerland)

For further information please contact:

Toni Meile, General Manager, Vascular International School AG, +41 79 611 01 04
meile@vascular-international.org