

SIMULATOR MODELS

Aorta | **Carotis** | Leg

Carotid simulation model for open carotid surgery

The VASCULAR INTERNATIONAL carotid simulator works independently and can be used anywhere. The water filled carotid inlay, under pulsatile pressure, makes a lifelike training of the conventional TEA and of the eversion-TEA possible.



Carotid simulation model ready to use



The new two-coloured carotid artery inlays work for interrupted sutures as well as continuous sutures.



Operation site after plaque removal



The whole carotid simulator is in a stackable simulator box which is easy to transport.



The carotid artery inlay is stenosed and makes a lifelike training of the conventional TEA and of the eversion-TEA possible.



All in one equipment

Each simulator comes with a full pump set including pump control and cable-tie kit. It can easily be transported in a box.

SALES PRICE incl. pump set, 1 carotid & 1 cable-tie kit CHF **2.950,-** (ex works Switzerland)

For further information please contact:

Toni Meile, General Manager, Vascular International School AG, +41 79 611 01 04
meile@vascular-international.org